Erba Lachema s.r.o.

ADVISORY NOTICE

EC 90 latest software 0.34.7 additional information

Cat. No.	Product
INS00079	EC 90

26th March 2021

Dear Valued Partner,

We would like to inform you that the current software that is installed on the EC 90 analyser has below listed deviations.
Please download and install the latest software version before performing electrolyte analysis.

The software version 0.34.7 currently installed on the analysers has the following known deviations:

HP measurement mode
If you select the High Precision mode for measurement, but do not aspirate the first sample package as requested the test results calculated will deviate from the real sample concentration.

Pop up window "measurement"
If you press any area besides the pop up window, after starting the measurement routine, the buttons of the surrounding area are still active and may cause unwanted process interactions (freeze).

In case you have any questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical Product Management

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -1-